
Miltex®

Redefining Excellence

# McGIVNEY HEMORRHOIDAL LIGATOR

## INSTRUCTIONS FOR USE


(Fig. 1)

The loading cone is inserted into the ligating drum of the instrument and the ligating o-ring is forced over the tip of the cone onto the surface of the instrument. (Fig. 1). Do not load the instrument until just before usage so that the elasticity of the ligating o-ring is not lost in prolonged stretching after being loaded on the instrument.

There are two thumbscrews on and near the handle and shaft of the instrument. The handle thumbscrew is used to assemble and secure the handle upon the shaft of the instrument. The thumbscrew which is either underneath or at the end of the shaft, depending on the model number, is used to facilitate rotation of the shaft in a 360-degree arc to make certain that the drum is in correct position


before ligating the hemorrhoid (Fig. 2). The thumbscrew on the shaft must be securely fastened after the proper rotation has been accomplished. Therefore the ligating drum may be placed in the proper quadrant of the lower rectum in order to facilitate the ligation of the hemorrhoid.

A fenestrated anoscope is inserted to expose the largest hemorrhoid which has been selected for treatment. The grasping forceps is held by the thumb and index finger in an inverted position with the hand in supination and the hemorrhoid pulled into the ligating drum, so that the o-ring may be placed around the neck of the hemorrhoid. The hemorrhoid is grasped sufficiently high enough to ensure the greatest amount of the hemorrhoid is grasped in the ligation. This point should be at least 10mm from the tip of the hemorrhoid towards the attachment


(Fig. 2)

(Fig. 3). If the hemorrhoid is grasped at its most dependent portion only a small amount of the hemorrhoid will be ligated and it will be at a point that is either on the anorectal line, or close enough to it to produce pain after the ligation. The ligation should be at least 7mm above the anorectal line. In this way the maximum amount of tissue is ligated. If the patient is asked to bear down at the moment of retracting the hemorrhoid and the instrument is forced against it, the amount of tissue ligated will be increased considerably (Fig. 4).

(Fig. 3)

When the maximum amount of hemorrhoid has been withdrawn into the ligating drum, the handle of the instrument is closed and the o-ring placed around the neck

(Fig. 4)

of the hemorrhoid. If the patient experiences pain at the time of the placement of the ligating o-ring, the hemorrhoid has been ligated too close to the anorectal line. The o-ring should be removed by means of a crypt hook and if necessary a pair of sharp, pointed scissors. A new o-ring is then applied at a higher level.

The patient may experience a minimum amount of pain shortly after the ligation but this soon disappears and most patients leave the office without any discomfort or inconvenience. Subsequently, a sense of fullness and a desire to defecate may be experienced by the patient, however, he should be cautioned to disregard this urge, otherwise the ligature o-ring could be passed at the time of the bowel movement and the procedure have to be repeated. This sensation disappears in a short time. The patient can be given a prescription for an analgesic to use in the event that there is sufficient discomfort to require relief. A written list of instructions can be provided to the patient to answer most of the questions he usually has after the procedure.

Usually the largest hemorrhoid present is ligated first and the smaller ones ligated at monthly intervals. No more than one hemorrhoid should be ligated at a time; otherwise, one increases the possibility of the complications of the procedure including thrombosis of the external hemorrhoidal plexus and adjacent hemorrhoids. A history of the patient should be taken to ascertain if the patient has a coagulation deficiency of the blood or is on any blood thinners before attempting treatment of the hemorrhoids. If this is done and the patient follows the instructions carefully, bleeding as a postoperative complication is negligible.

Following the ligation, the hemorrhoid swells, becomes dark in color and within a short time infarction occurs followed by necrosis and sloughing of the hemorrhoid. Bleeding following the sloughing of the hemorrhoid may occur in small amounts. Should this become excessive, the patient should notify the doctor immediately so that measures may be taken to stop the hemorrhage. The bleeding may be controlled by means of cryotherapy or electrosurgery in the office or the outpatient clinic. Since the anal canal has not been invaded by the operation, anoscopy is painless and one may examine the patient without discomfort and with full cooperation. Those patients who are accustomed to taking laxatives should continue to do so in order to assure a daily bowel movement to avoid fecal impaction. Others are cautioned to take special care to have a daily bowel movement and to avoid enemas.

The instrument should be cleaned and sterilized after each use and periodically disassembled and lubricated with a non-silicone lubricant such as Miltex Spray Lube.

*Latex O-rings (#28-151) contain natural rubber latex which may cause allergic reactions in some individuals.*

*Latex-free O-rings (#28-155) are also available.*

### THE McGIVNEY HEMORRHOIDAL LIGATOR - Designed By John McGivney, M.D., F.A.C.S.

REF	Description
28-150A	McGIVNEY Hemorrhoidal Ligator Set - includes McGIVNEY Hemorrhoidal Ligator 7" (17.8 cm) shaft, McGIVNEY Grasping Forceps, Loading Cone and Ligating O-Rings. Compatible with LATEX O-Rings (PN 28-151) only.
28-150	McGIVNEY Hemorrhoidal Ligator, 7" (17.8 cm) shaft. Includes Loading Cone and Ligating O-Rings. Compatible with LATEX O-Rings (PN 28-151) only.
28-150XL	McGIVNEY Hemorrhoidal Ligator, 10" (25.4 cm) shaft. Includes Loading Cone and Ligating O-Rings. Compatible with LATEX O-Rings (PN 28-151) only.
28-151	LATEX Ligating O-Rings compatible with PN's 28-150, 28-150XL & 28-154B Ligator patterns. 100 rings per bag. Provided Non-Sterile.
28-152	McGIVNEY Grasping Forceps, 7-1/2" (19.1 cm)
28-153	Loading Cone compatible with PN 28-150 & 28-150XL Ligators only.
28-154A	McGIVNEY Hemorrhoidal Ligator Set for LATEX-FREE applications - includes McGivney Hemorrhoidal Ligator 8" (20.3 cm) shaft, McGIVNEY Grasping Forceps, Loading Cone and LATEX-FREE Ligating O-Rings. Compatible with LATEX-FREE & Standard LATEX O-Ring PN's 28-151 & 28-155 only.
28-154B	McGIVNEY Hemorrhoidal Ligator for LATEX-FREE applications, 8" (20.3 cm) shaft. Includes Loading Cone and LATEX-FREE Ligating O-Rings. Compatible with LATEX-FREE & Standard LATEX O-Ring PN's 28-151 & 28-155 only.
28-155	LATEX-FREE Ligating O-Rings compatible with PN 28-154B Ligator only. 100 rings per bag. Provided Non-Sterile.
28-156	Loading Cone compatible with PN 28-154B Ligator only.

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