


# MILTEX HI-LIGHT


## MICRO MIRRORS



**Rhodium Front Surface Mirror**  
Provides "true to life" images.

**Small and Flexible Mirror Head**  
Allows visualization of hard to see areas.

**Variety**  
Available in four different sizes and shapes.



The Miltex® Micro Mirrors feature a rhodium front surface mirror fitted on a stainless steel mounting to assure sharpness and clarity of the reflected image. Other micro mirrors on the market with highly polished stainless steel surfaces are incomparable to the Miltex® Micro Mirrors in producing high quality images. The mirrors are used in conjunction with microscopes during micro surgery such as periodontitis and endodontic surgeries. The compact size and flexibility of the mirror head allow a very close approach to the most inaccessible points under examination.

REF	Description
017-24800	MICRO MIR 3MM OVAL CS FS
017-24801	MICRO MIR 5MM OVAL CS FS
017-24802	MICRO MIR 9/3MM REC CS FS
017-24803	MICRO MIR 7/2MM REC CS FS

*Miltex*  
Redefining Excellence

phone 717 840-9335  
toll free 866 854-8300  
www.miltex.com

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.