


# MILTEX® HI-LIGHT

## MILTEX® MICROSURGICAL INSTRUMENTS



**REF 6890** \*with MILTEX® Storage Cassette & Mat

**MILTEX® Microsurgical Kit Includes 1 each of:**

8-4	DERF Needle Holder	grated jaws, 4000 teeth per square inch	4-3/4" (12.1cm)
5-300	Iris Scissors	straight	4" (10.2 cm)
17-305X	SWISS Jeweler Style Forceps	non-magnetic stainless steel, style 5, micro-fine	4-3/8" (11.2 cm)
18-948	McPHERSON Micro Iris Suturing Forceps	straight, smooth platform 0.3 mm X 5 mm, no teeth	3-1/2" (8.9 cm)
18-1625	McPHERSON-VANNAS Micro Iris Scissors	straight, sharp points	3-1/4" (8.3 cm)
18-34	Plain Wire Eye Speculum	blades 15 mm wide	2-1/2" (6.4 cm)
10-2-ST	MILTEX® Probe	double ended, malleable, sterling	4-1/2" (11.4 cm)
4-6835	MILTEX® Storage Cassette	dimensions 6-1/2" x 6-1/4" x 1-1/4"	-
4-6835MAT	MILTEX® Storage Cassette Mat	-	-

**REF 6890-NC** \*without MILTEX® Storage Cassette & Mat

See reverse side for more microsurgical instruments

Illustrations and content provide general description only and may be subject to change. For additional information, to schedule an in-office demonstration, or to place an order, please contact your local Miltex authorized Distributor or Miltex Customer Service at 866-854-8300, e-mail: [customerservice@miltex.com](mailto:customerservice@miltex.com) or visit us at [www.miltex.com](http://www.miltex.com).

Miltex is a registered trademark of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation.  
The Integra wave logo is a trademark of Integra Life Sciences Corporation. Copyright © 2008 Integra LifeSciences Corporation

*Miltex*  
Redefining Excellence

**INTEGRA**

toll free 866 854-8300  
fax 866 854-8400  
[www.miltex.com](http://www.miltex.com)



## Scissors

REF	Description	Type	Size
17-2060	MILTEX® Microsurgical Scissors	sharp points, straight, 8 mm blades, flat handles	6" (15.2 cm)
17-2040	MILTEX® Microsurgical Scissors	sharp points, curved, 8 mm blades, flat handles	6" (15.2 cm)
17-2150	MILTEX® Microsurgical Scissors	sharp points, straight, 6 mm blades, round handles	6" (15.2 cm)
17-2100	MILTEX® Microsurgical Scissors	sharp points, curved, 6 mm blades, round handles	6" (15.2 cm)
18-1618	Micro Iris Scissors	straight, sharp points	4" (10.2 cm)
18-1619	Micro Iris Scissors	curved, sharp points	4" (10.2 cm)
18-1625	MCPHERSON-VANNAS Micro Iris Scissors	straight, sharp points	3-1/4" (8.3 cm)
18-1626	MCPHERSON-VANNAS Scissors	curved, sharp points	3-1/4" (8.3 cm)
18-1630	VANNAS Micro Scissors	straight, ultra fine blades	3-1/4" (8.3 cm)
18-1631	VANNAS Micro Scissors	curved, ultra fine blade	3-1/4" (8.3 cm)

17-2060



## Forceps

REF	Description	Type	Size
17-301	SWISS Jeweler Style Forceps	non-magnetic stainless steel, style 1, fine	4-3/4" (12.1 cm)
17-302	SWISS Jeweler Style Forceps	non-magnetic stainless steel, style 2, strong, fine	4-3/4" (12.1 cm)
17-303	SWISS Jeweler Style Forceps	non-magnetic stainless steel, style 3, narrow, fine	4-3/4" (12.1 cm)
17-303C	SWISS Jeweler Style Forceps	non-magnetic stainless steel, style 3C, narrow, fine	4-3/8" (11.2 cm)
17-303X	SWISS Jeweler Style Forceps	non-magnetic stainless steel, style 3, micro-fine	4-3/4" (12.1 cm)
17-304	SWISS Jeweler Style Forceps	non-magnetic stainless steel, style 4, fine	4-3/8" (11.2 cm)
17-305	SWISS Jeweler Style Forceps	non-magnetic stainless steel, style 5, super fine	4-3/8" (11.2 cm)
17-305X	SWISS Jeweler Style Forceps	non-magnetic stainless steel, style 5, micro-fine	4-3/8" (11.2 cm)
17-306	SWISS Jeweler Style Forceps	non-magnetic stainless steel, style 6, pickup	4-1/2" (11.4 cm)
17-307	SWISS Jeweler Style Forceps	non-magnetic stainless steel, style 7, fine, curved	4-1/2" (11.4 cm)

17-301



## Needle Holders

REF	Description	Type	Size
8-4TC	DERF Needle Holder	grated jaws, 4000 teeth .per square inch	4-3/4" (12.1cm)
8-14TC	OLSEN-HEGAR Needle Holder with Suture Scissors	extra delicate, serrated jaws, 4000 teeth per square inch	4-3/4" (12.1 cm)
17-1000	MILTEX® Microsurgical Needle Holders	round handles, 0.6 mm tips, curved jaws	5-1/4" (13.3 cm)
17-1010	MILTEX® Microsurgical Needle Holders	round handles, 0.6 mm tips, curved jaws	7-1/8" (18.1 cm)
17-1020	MILTEX® Microsurgical Needle Holders	round handles, 0.6 mm tips, curved jaws, with lock	5-1/4" (13.3 cm)
17-1024	MILTEX® Microsurgical Needle Holders	round handles, 0.6 mm tips, straight jaws, with lock	7-1/8" (18.1 cm)
17-1025	MILTEX® Microsurgical Needle Holders	round handles, 0.6 mm tips, curved jaws, with lock	7-1/8" (18.1 cm)
17-1090TC	JACOBSON Micro Vascular Needle Holder	straight, with lock, Carb-N-Sert jaws, 16000 teeth P. Sq. inch	7" (17.8 cm)
17-1091TC	JACOBSON Micro Vascular Needle Holder	curved, with lock, Carb-N-Sert jaws, 16000 teeth P. Sq. inch	7" (17.8 cm)
17-1092TC	JACOBSON Micro Vascular Needle Holder	straight, with lock, Carb-N-Sert jaws, 16000 teeth P. Sq. inch	8-1/2" (21.6 cm)

18-4TC



## Retractors

REF	Description	Type	Size
18-2	LANCASTER Eye Speculum	blades 15 x 5 mm	3-1/4" (8.3 cm)
18-4	KNAPP Eye Speculum	blades 15 x 6 mm	3-1/4" (8.3 cm)
18-10	CASTROVIEJO Eye Speculum	large blades 18 x 5 mm	4" (10.2 cm)
18-22	WILLIAM Eye Speculum	large blades, 15 x 5 mm	3-1/4" (8.3 cm)
18-24	WILLIAMS Eye Speculum	small blades, 11 x 5 mm	2-3/4" (7 cm)
18-35	BARRAQUER Eye Speculum	large blades 15 mm wide	1-3/4" (4.4 cm)
18-46	TEMPORAL APPROACH Wire Speculum	spring on nasal side, solid blades	1-1/2" (3.8 cm)
18-47	SAUER Eye Speculum	petite size, blades 10 mm wide	1-1/4" (3.2 cm)
18-50	MURDOCK Eye Speculum	blades 15 mm wide	2-1/4" (5.7 cm)
18-54	WIENER Eye Speculum	with extensions to cover eye lashes, blades 15 mm	2-1/4" (5.7cm)
18-63	COOK Eye Speculum	petite size, 9 mm blades	1-7/8" (4.7 cm)
18-66	MELLINGER Eye Speculum	solid blades 15 mm blades	2-3/4" (7 cm)
18-80	AGRICOLA Lacrimal Sac Retractor	3 x 3 pointed prongs 4 mm deep	1-1/2" (3.8 cm)
18-84	MELLER Lacrimal Sac Retractor	3 x 3 prongs 3 mm deep	2-7/8" (7.3 cm)
18-88	GOLDSTEIN Lacrimal Sac Retractor	3 x 3 point prongs 5 mm deep	1-1/8" (2.8 cm)

18-2



## Punch


REF	Description
76-5	AINSWORTH Punch

76-5



## Biopsy Punch with Plunger System

REF	Size
33-31AA-P/25	1.0 mm w/plunger, sterile, disposable
33-31A-P/25	1.5 mm w/plunger, sterile, disposable
33-31-P/25	2.0 mm w/plunger, sterile, disposable



33-31AA-P/25

33-31A-P/25

33-31-P/25