

MILTEX HI-LIGHT

SPRINT TCM



Confidently clean and shape curved canals with the Sprint Torque Controlled Motor. The Sprint TCM has preselectable speed and torque limits, greatly reducing the risk of file separation. When the NiTi file encounters the preset torque limit level, the Sprint TCM will stop and reverse direction. Step up to the next generation in endodontic motors with the Sprint TCM.

FEATURES

- Two selectable torque modes, greatly reducing the chance of file separation:
 1. Auto stop
 2. Auto reverse
- Full power and selectable torque levels from 2 to 50gm/cm
- Three reduction angles (1:1, 8:1 and 10:1)
- Full low-speed range 1200 to 16000 rpms for consistent speed control

Sets:

REF	Description
017-36295	Sprint TCM and Tardie™ 10:1 Handpiece
017-36296	Sprint TCM and Tardie™ 8:1 Handpiece

Components:

REF	Description
017-36300	Sprint TCM
017-36307	Tardie™ 10:1 Reduction Handpiece
017-36308	Tardie™ 8:1 Reduction Handpiece
017-36309	Tardie™ 1:1 Handpiece

Miltex
Redefining Excellence

589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.

0704/5/B/1

SPRINTTCM-HI

