


# Integra™ Miltex®

Veterinary Surgical and Dental  
Instrument Packs

Limit uncertainty with  
convenience and  
ease of use



**INTEGRA**   
LIMIT UNCERTAINTY

**Limit uncertainty with Integra™ Miltex® surgical and dental instrument kits which have been developed through the guidance of veterinary professionals.**

**Our veterinary instrument packs are designed for convenience, ease of use and are organized in stainless steel sterilization and storage cassettes.**

We offer kits in the following instrument grades:

**Miltex® Premium Grade:**

- German stainless steel
- German forgings
- German crafted finishing

**Miltex® Standard Grade:**

- German stainless steel
- German forgings
- European finishing

**Miltex® Mid Grade:**

- Japanese stainless steel
- Pakistani forgings
- Pakistani finishing
- Cost competitive

**Request our comprehensive Animal Health Catalog by contacting our Customer Service Team at 866-854-8300 or [customerservice.miltex@integralife.com](mailto:customerservice.miltex@integralife.com)**




### 6800 - General Surgery Kit

Kit includes:

Qty.	Reference	Description
1	5-124	Mayo Scissors, Straight, 6-3/4" (17.1 cm)
1	4-7	#3 Scalpel Handle
1	6008	Snook Hook, 8" (20.3 cm)
3	7-4	Mosquito Forceps, Curved, 5" (12.7 cm)
2	7-36	Kelly Forceps, Straight, 5-1/2" (14 cm)
1	5-182	Metzenbaum Scissors, Curved, 7" (17.8 cm)
1	6-124	Brown-Adson Forceps, 7 X 7 Teeth, 4-3/4" (12.1 cm)
4	7-172	Rochester-Carmalt Forceps, Curved, 6-1/4" (15.9 cm)
4	7-506	Backhaus Towel Clamps, 5-1/4" (13.3 cm)
1	8-15TC	Olsen-Hegar Needle Holder, Tungsten Carbide, 5-1/2" (14 cm)
1	3-6815	Storage Cassette, 15" (38.1 cm) X 8-1/4" (21 cm) X 1-1/4" (3.1 cm)

	Reference
<b>Premium Grade</b>	6800
<b>Standard Grade</b>	MH6800
<b>Mid Grade</b>	V96800



### 6810 Canine / 6820 Feline - Spay Kit

Kit includes 1 each of:

Qty.	Reference	Description
1	4-7	#3 Scalpel Handle
1	5-124	Mayo Scissors, Straight, 6-3/4" (17.1 cm)
1	6008	Snook Hook, 8" (20.3 cm)
2	7-38	Kelly Forceps, Curved, 5-1/2" (14 cm)
1	5-182	Metzenbaum Scissors, Curved, 7" (17.8 cm)
1	6-44	Tissue Forceps, 1 X 2 Teeth, 5-1/2" (14 cm)
1	16-6	Allis Tissue Forceps, 4 X 5 Teeth, 6" (15.2 cm)
4	7-504	Backhaus Towel Clamps, 3-1/2" (8.9 cm)
1	6-124	Brown-Adson Forceps, 7 X 7 Teeth, 4-3/4" (12.1 cm)
1	8-16TC*	Olsen-Hegar Needle Holder, Tungsten Carbide, 6-1/2" (16.5 cm)
	or	
1	8-15TC*	Olsen-Hegar Needle Holder, Tungsten Carbide, 5-1/2" (14 cm)
2	7-4	Halsted Mosquito Forceps, Curved, 5" (12.7 cm)
2	7-172	Rochester-Carmalt Forceps, Curved, 6-1/4" (15.9 cm)
1	3-6815	Storage Cassette, 15" (38.1 cm) X 8-1/4" (21 cm) X 1-1/4" (3.1 cm)

	Reference - Canine	Reference - Feline
<b>Premium Grade</b>	6810	6820
<b>Standard Grade</b>	MH6810	MH6820
<b>Mid Grade</b>	V96810	V96820




\* The canine kit includes an 8-16TC; feline kit includes an 8-15TC

### 6830 - Basic Veterinary Eye Kit

Kit includes:

Qty.	Reference	Description
1	4-7	#3 Scalpel Handle
1	18-35	Wire Eye Speculum
1	18-780	Eye Dressing Forceps, 4" (10.2 cm)
1	18-1462	Stevens Tenotomy Scissors, Straight, 4-1/8" (10.5 cm)
1	18-1202	Chalazion Forceps, 3-1/2" (8.9 cm)
1	18-1828	Castroviejo Needle Holder with Lock, Straight, 5-1/2" (14 cm)
1	3-6835	Storage Cassette, 6-1/2" (16.5 cm) X 6-1/4" (15.8 cm) X 1-1/4" (3.1 cm)



**6840 - Orthopedic Kit**

Kit Includes:

Qty.	Reference	Description
1	9-40-20	Babcock Stainless Steel Suture Wire (Gauge # 20)
1	9-40-30	Babcock Stainless Steel Suture Wire (Gauge # 30)
1	9-124	Wire Cutting Scissors, One Serrated Blade, 4-3/4" (12.1 cm)
1	11-74	Senn Retractor, Sharp, Double Ended, 6-3/8" (16.2 cm)
1	11-76	Senn Retractor, Blunt, Double Ended, 6-3/8" (16.2 cm)
1	11-608	Weitlaner Retractor, Sharp, 3 X 4 Teeth, 5-1/2" (14 cm)
1	18-1968	Freer Elevator, 7" (17.8 cm)
1	19-708	Spratt Curette, Size 1/0 (3.7 X 5.8 mm), 6-3/4" (17.1 cm)
1	19-792	Lucae Mallet, Solid Stainless, 7-3/4" (19.7 cm)
1	21-369	Fomon Rasp, Double Ended, Coarse, 8-1/4" (21 cm)
1	22-500TC	Wire Twister, Tungsten Carbide, 6" (15.2 cm)
2	26-120	Gigli-Strully Saw Handle, Standard Pattern, T-Shaped
1	26-130	Gigli Saw, Standard Twisted Wire, 12" (30.5 cm)
1	27-4	Kern Bone Holding Forceps, with Ratchet, 9-1/2" (24.1 cm)
1	27-7	Kern Bone Holding Forceps, with Ratchet, 6-3/4" (17.1 cm)
1	27-158	Cannulated Pin and Wire Cutter, 7-1/2" (19.1 cm)
1	27-1390	Liston Bone Cutting Forceps, Straight, 5-1/2" (14 cm)
1	27-603	Swiss Pattern Osteotome, 14 mm Wide
1	3-6845	Storage Cassette, 10-3/4" (27.30 cm) X 8-1/4" (21 cm) X 2-3/4" (7 cm)



**27-610 - Osteotome Kit**

Kit includes:

**Description**  
 Set of 9 Swiss Pattern Osteotomes, 2 mm to 20 mm, Complete in Customized Storage Cassette (10-1/2" (26.6 cm) X 8" (20.3 cm) X 1/2" (1.2 cm))

- 2 mm
- 4 mm
- 6 mm
- 8 mm
- 10 mm
- 12 mm
- 14 mm
- 18 mm
- 20 mm



### 6890 - Microsurgical Instrument Kit (with Cassette and Mat)

### 6890-NC - Microsurgical Instrument Kit (without Cassette and Mat)

Kit Includes:

Qty.	Reference	Description
1	8-4	Derf Needle Holder, Serrated Jaws, 4-3/4" (12.1 cm)
1	5-300	Iris Scissors, Straight, 4" (10.2 cm)
1	17-305X	Swiss Jeweler Style Forceps, Style 5, Micro-Fine, 4-3/8" (11.2 cm)
1	18-948	McPherson Micro Iris Suturing Forceps, Straight, 3-1/2" (8.9 cm)
1	18-1625	McPherson-Vannas Micro Iris Scissors, Straight, Sharp Points, 3-1/4" (8.3 cm)
1	18-34	Plain Wire Speculum, Blades 15 mm Wide, 2-1/2" (6.4 cm)
1	10-2-ST	Probe, Double Ended, Malleable, Sterling, 4-1/2" (11.4 cm)
1	4-6835	Storage Cassette, 6-1/2" (16.5 cm) X 6-1/4" (15.8 cm) X 1-1/4" (3.1 cm)
1	4-6835MAT	Mat For Storage Cassette 4-6835




### 6891 - Lab Animal Kit (with Storage Cassette)

### 6891-NC - Lab Animal Kit (without Storage Cassette)

Kit includes:

Qty.	Reference	Description
1	5-181	Metzenbaum Scissors, Straight, 7" (17.7 cm)
1	5-284	Metzenbaum Scissors, Curved, 5" (12.7 cm)
1	6-42	Tissue Forceps, 1 X 2 Teeth, Serrated Handles, 5" (12.7 cm)
1	6-6	Dressing Forceps, Serrated Handles, 5" (12.7 cm)
1	18-956	Nugent Utility Forceps, 45° Angled Jaws, 4-1/4" (10.7 cm)
1	22-481	Blumenthal Oral Rongeur, Jaws at 30° Angle, 7" (17.7 cm)
1	4-6835	Storage Cassette, 6-1/2" (16.5 cm) X 6-1/4" (15.8 cm) X 1-1/4" (3.1 cm)



**6809 - Weldin Periodontal Instrument Kit**

Kit includes:

Qty.	Reference	Description
1	SVET-5450	Griplite® S6 Mirror Handle
1	67-696/3*	#3 Cone Socket Front Surface Mirror
1	67-696/4*	#4 Cone Socket Front Surface Mirror
1	SVET-23E/18P	Griplite® S6 #23 Explorer Probe
1	SVET-5004	Griplite® S6 Jacquette J33/1Ys Scaler
1	SVET-767	Griplite® S6 Morse Scaler
1	SVET-5201	Griplite® S6 Columbia, 13/14 Universal Curette
1	SVET-5202	Griplite® S6 Barnhart, 1/2 Universal Curette
1	DTST**	Delron Sharpening Test Stick
1	STN-ARK-CON	Conical Arkansas Sharpening Stone
1	STN-ARK-FLT-4	#4 Flat Arkansas Sharpening Stone
1	3-084107	Storage Cassette, 8" (20.3 cm) X 4-7/8" (10.1 cm) X 1-1/8" (2.5 cm)
1	-	Training DVD & Instructional Guide

\*(not sold individually - reference is 12/bx)      \*\*(not sold individually - reference is 6/pk)



**70-177 - Periotome Kit**

Ideal for separating the periodontal ligament fibers from the tooth in preparation for a dental extraction.

Kit Includes:

Qty.	Reference	Description
1	70-170	Periotome, Straight, Double Ended
1	70-171	Periotome, Hatchet, Double Ended
1	70-172	Periotome, Hoe, Double Ended
1	3-072007	Storage Cassette, 7" (17.7 cm) X 2-5/8" (5 cm) X 11/16" (1.9 cm)



**6806 - Extraction Instrument Kit**

Kit includes:

Qty.	Reference	Descriptions
1	4-7	# 3 Scalpel Handle
1	73-64	Baby Molt Periosteal Elevator, 6" (15.2 cm)
1	DEFF	Extracting Forceps Universal Mini, 4-1/2" (11.4 cm)
1	5D-324	La Grange Scissors, Curved, One Serrated Blade, Angled Shank, 4-1/4" (10.8 cm)
1	8-15TC	Olsen-Hegar Needle Holder, Tungsten Carbide, 5-1/2" (14 cm)
1	6-124	Brown-Adson Forceps, 7 X 7 Teeth, 4-3/4" (12.1 cm)
1	5-16	Operating Scissors, Straight, Sharp/Blunt, 5-1/2" (14 cm)
1	DELLUXS2	Luxating Elevator, Straight, 2 mm, 4-1/2" (11.4 cm)
1	DELLUXS3	Luxating Elevator, Straight, 3 mm, 4-1/2" (11.4 cm)
1	DELLUXS5	Luxating Elevator, Straight, 5 mm, 4-1/2" (11.4 cm)
1	DELLUXC3	Luxating Elevator, Curved, 3 mm, 4-1/2" (11.4 cm)
1	DELLUXC5	Luxating Elevator, Curved, 5 mm, 4-1/2" (11.4 cm)
1	3-085009	Storage Cassette, 8" (20.3 cm) X 6" (15.2 cm) X 1-3/16" (2.5 cm)



### DELLUXKIT - Luxating Dental Elevator Kit

The Integra™ Miltex® Luxating Dental Elevator Kit is a complete set of elevators ranging in size from 2 mm through 5 mm with curved and straight tip designs. The Luxators are designed to cut periodontal ligament while easing the tooth from the socket by rocking and cutting, rather than by lifting and prying like a conventional elevator.

**Kit Includes:**

Qty.	Reference	Description
1	DELLUXC2	Luxating Elevator, Curved, 2 mm, 4-1/2" (11.4 cm)
1	DELLUXC3	Luxating Elevator, Curved, 3 mm, 4-1/2" (11.4 cm)
1	DELLUXC4	Luxating Elevator, Curved, 4 mm, 4-1/2" (11.4 cm)
1	DELLUXC5	Luxating Elevator, Curved, 5 mm, 4-1/2" (11.4 cm)
1	DELLUXS2	Luxating Elevator, Straight, 2 mm, 4-1/2" (11.4 cm)
1	DELLUXS3	Luxating Elevator, Straight, 3 mm, 4-1/2" (11.4 cm)
1	DELLUXS4	Luxating Elevator, Straight, 4 mm, 4-1/2" (11.4 cm)
1	DELLUXS5	Luxating Elevator, Straight, 5 mm, 4-1/2" (11.4 cm)
1	3-084007	#7 Storage Cassette, 8" (20.3 cm) X 4-7/8" (10.1 cm) X 1-1/8" (2.5 cm)



### DELWDELUX - Winged Dental Elevator Deluxe Kit

**Kit includes:**

Qty.	Reference	Description
1	DELWNG2	Winged Elevator, 2 mm, 4-1/2" (11.4 cm)
1	DELWNG3	Winged Elevator, 3 mm, 4-1/2" (11.4 cm)
1	DELWNG4	Winged Elevator, 4 mm, 4-1/2" (11.4 cm)
1	DELWNG5	Winged Elevator, 5 mm, 4-1/2" (11.4 cm)
1	DELWNG6	Winged Elevator, 6 mm, 4-1/2" (11.4 cm)
1	DELWNG7	Winged Elevator, 7 mm, 4-1/2" (11.4 cm)
1	DELWNG8	Winged Elevator, 8 mm, 4-1/2" (11.4 cm)
1	3-084007	#7 Storage Cassette, 8" (20.3 cm) X 4-7/8" (10.1 cm) X 1-1/8" (2.5 cm)



Winged Dental Elevators are shaped to improve instrument adaptation to the root structure. When used in a side-to-side rocking/twisting motion, the blades provide additional leverage, especially when used between teeth in the interdental space. The winged-blades will hook under the buccal bulge, improving the angle of direction in order to lift the tooth and root.

### DELWBASIC - Winged Dental Elevator Basic Kit

**Kit includes:**

Qty.	Reference	Description
1	DELWNG2	Winged Elevator, 2 mm, 4-1/2" (11.4 cm)
1	DELWNG3	Winged Elevator, 3 mm, 4-1/2" (11.4 cm)
1	DELWNG4	Winged Elevator, 4 mm, 4-1/2" (11.4 cm)
1	DELWNG5	Winged Elevator, 5 mm, 4-1/2" (11.4 cm)
1	3-083005	#5 Storage Cassette, 8" (20.3 cm) X 3-13/16" (7.6 cm) X 1-1/8" (2.5 cm)





\*+H834VETPACKHIOY\*

For more information or to place an order, please contact:  
Integra York PA, Inc. ■ 589 Davies Drive, York, PA 17402  
866-854-8300 USA ■ 717-840-2763 outside USA ■ 717-840-9335 fax  
[integralife.com/integra-miltex](http://integralife.com/integra-miltex)

