


MILTEX HI-LIGHT

PERIOTOMES

AVAILABLE IN THREE DOUBLE-ENDED STYLES
STRAIGHT-BLADE, HATCHET-BLADE AND HOE-BLADE



Miltex® Periotomes are designed for separating the periodontal ligament fibers from the tooth in preparation for a dental implant. After the periodontal ligament is cut and the tooth removed, the underlying bone structure remains intact, giving the necessary support for an implanted replacement tooth.

The Periotome Kit contains one of each pattern in a Slimline Cassette and is an excellent way to store, organize, and protect Periotomes.

PERIOTOME KIT (70-177)

CONTAINS ONE OF EACH:

70-170	Straight Blade
70-171	Hatchet Blade
70-172	Hoe Blade
3-082008	Slimline Cassette 8" x 2-5/8" x 11/16"



70-170
Straight Blade



70-171
Hatchet Blade



70-172
Hoe Blade

Miltex
Redefining Excellence

589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.