


MILTEX HI-LIGHT

JUNIORS DISPOSABLE ENDODONTIC MEDICAMENTS



Precise dispensing every time with the new Miltex® disposable endodontic medicaments. The squeezable bulb eliminates the need to transfer solutions from container to syringe.

FEATURES

- Color-coded pipettes allow for easy identification of solutions
- Squeezable bulbs offer accuracy similar to syringes
- Angled stems allow easy access into canals

REF	Product Name	Description
033-39210P	Slide Juniors	File lubricant w/10% urea peroxide
033-39211P	Slide SH Juniors	Slide w/3% Sodium Hypochlorite
017-39211P	Endodilator Juniors	15% EDTA
033-39212P	Chorhexidine Juniors	2% Chlorhexidine
033-39214P	Endo-Klean Juniors	2.5% Sodium Hypochlorite
033-39215P	Lube Juniors	File lubricant w/NaOH and EDTA

Miltex
Redefining Excellence

phone 717 840-9335
toll free 866 854-8300
www.miltex.com

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.

