


MILTEX® HI-LIGHT

MEMBRANE TACK KIT



62-305
(3 mm)

62-306
(5 mm)

62-300

The Miltex® Membrane Tack Kit provides dental clinicians with the proper instrumentation for fixating or stabilizing guided tissue regeneration membranes to bone. The Membrane Tacks are small sharp titanium alloy pins used during the healing process by providing an attachment mechanism for the membrane to local and adjacent bone at the surgical site.

The membrane tack kit is easy to use. The membrane tack applicator reduces handling from storage box to placement against membrane in one easy step. (For detailed placement instructions, please refer to Instructions for Use included with the product.)

Note: Membrane tacks are not intended for long-term implantation, but should be removed upon completion of the healing process.

Benefits of Membrane Tacks

- High Quality German Titanium membrane tacks are strong and composed of corrosion-resistant metal
- Available in two sizes: (5mm, 3mm)
- Easy delivery to surgical site by simply engaging tack applicator onto membrane tack for hands free placement

Benefits of Membrane Tack Kit

- High quality German stainless steel instruments and cassette
- Membrane Tack Kit cassette stores, organizes and protects all instruments; autoclavable
- Convenient Membrane Tack storage box holds up to 15 tacks, providing enough tacks per surgical case

See reverse side for Membrane Tack Kit details

MILTEX® MEMBRANE TACK KIT

The membrane tack kit contains the following:

REF	Description
62-300	Dental Membrane Tack Kit
62-301	Membrane Tack Probe
62-302	Membrane Tack Applicator
62-303	Perforation Raspatory
62-304	KIRSCH, Sinus 7, Double-ended
62-305	Membrane Tack - Titanium 3mm 5/Pk
62-306	Membrane Tack - Titanium 5mm 5/Pk
62-307	Sinus Elevator 2
62-308	Mixing Cup, Stainless Steel
62-309	Membrane Tack Storage Box
62-311	Kit Tack Wash Tray with Lid

MEMBRANE TACK KIT ACCESSORIES

REF	Description
62-310	Mallet 16.5 cm



62-301 Membrane Tack Probe



62-302 Membrane Tack Applicator



62-303 Perforation Raspatory



62-304 Double-ended Sinus 7 Instrument



62-305 Membrane Tack - Titanium 3mm 5/Pk



62-306 Membrane Tack - Titanium 5mm 5/Pk



62-307 Sinus Elevator



62-308 Mixing Cup with Lid



62-309 Membrane Tack Storage Box
(Membrane tacks not included)



62-311 Kit Wash Tray with Lid
(Contents not included)



62-310 Mallet 16.5 cm



For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.

Illustrations and content provide general description only and may be subject to change. Miltex and Miltex Redefining Excellence are registered trademarks of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation. The Integra wave logo is a trademark of Integra LifeSciences Corporation. Copyright © 2010 Integra LifeSciences Corporation