

Integra™ Miltex®

Animal Health

BIANNUAL SPECIALS
January 1 – June 30, 2012



INTEGRA 
LIMIT UNCERTAINTY

Integra™ Miltex® Loupes and LED Light System

Limit uncertainty with portability and amazing brightness.

- Integra™ Miltex® optics quickly and easily lock into your personal setting.
- Get a custom fit without the need for individual measurement!
- LED light has an innovative power pack design with external battery cells for quick changing.
- High intensity portable dental and surgical headlight (80,000 Lux @ 14 inches).
- Affixes easily to your Integra™ Miltex® Loupes.



A New CryoSurgery Solution!

Purchase a complete set (33510) and receive an additional tip at no charge!

CryoSolutions® is a cryosurgery device that uses liquefied Nitrous Oxide at an -89°C (-127°F) working temperature. Cryosurgery has been used in veterinary dermatology for several decades and is used to treat many types of skin lesions. The advantages of cryosurgery include easy application, short treatment time, low pain, high success rate, infrequent skin damage, infrequent scarring and modest treatment costs. This product provides a consistent and ideal temperature for localized application to the skin lesion without damaging surrounding tissue.

CryoSolutions, Product Listing and Replacement Parts

Reference	Description
33510	Complete set, includes unit with standard 1 mm wide tip, one cartridge (23.5 g N2O), metal pin, user manual and plastic case
33511	Standard 1 mm wide glass tip
33512	Dermatology 2 mm wide glass tip
33513	Dermatology 3 mm wide glass tip
33514	Dermatology 4 mm wide glass tip
33516	10 Pack of Cartridges (23.5 g N2O) medical grade, with threaded safety valve
33517	4 Pack of Cartridges (23.5 g N2O) medical grade, with threaded safety valve
33519	CryoSolutions replacement plastic carrying case with protective foam
33520	CryoSolutions replacement metal pin



Fluoride Coated Nurse, Bandage and Utility Scissors Buy 3, get 1 Free (mix and match)*

Nurse Scissors, Fluoride Coated with colored safety guard

Reference	Description
5703	Nurse Scissors 5-1/2" Black
5704	Nurse Scissors 5-1/2" Neon Green
5705	Nurse Scissors 5-1/2" Lavender



Bandage and Utility Scissors, Fluoride Coated with serrated lower blade

Reference	Description
5803	Bandage & Utility Scissors 8" Hot Pink
5804	Bandage & Utility Scissors 8" Bright Blue
5801	Bandage & Utility Scissors 8" Black
5802	Bandage & Utility Scissors 8" Lavender





* All buy/get promotions apply only to Integra™ Miltex® products advertised in this quarterly flyer unless otherwise noted.

* Redeemed goods must be of equal or lesser value of lowest priced purchased item.

Curettes

Buy 3, get 1 (mix and match)*

Reference	Description
80-5203	Barnhart 5/6
80-5204	McCalls 13S/14S
SVET-5201	Columbia 13/14
SVET-5202	Barnhart 1/2




SVET-5202

Scalers

Buy 3, get 1 (mix and match)*


Reference	Description
SVET-767	Morse O-00
SVET-5004	Jacquette 33/H5



Probes/Explorers

Buy 3, get 1 (mix and match)*

Reference	Description	Markings
SVET-23E/18P	18 Probe - 23	1, 2, 3, 5, 6, 7, 8, 10, 11, 12, 13, 15, 16, 17, 18 mm
80-5401	OW Probe - 23 Explorer	1, 2, 3, 5, 7, 8, 9, 10 mm



Buy \$250 in dental instrumentation and receive:

- (1) #4 Flat Arkansas Conical Sharpening Stone (STN-ARK-FLT-4)
- (1) Arkansas Conical Sharpening Stone (STN-ARK-CON)
- (1) Sharpening Oil (STN-OIL)



* All buy/get promotions apply only to Integra™ Miltex® products advertised in this quarterly flyer unless otherwise noted.

* Redeemed goods must be of equal or lesser value of lowest priced purchased item.

Dental Kits

Buy 3, get 1 (mix and match)*

6809 - Weldin Periodontal Instrument Kit

Kit includes:

Qty.	Reference	Description
1	SVET-5450	Griplite® S6 Mirror Handle
1	67-696/3*	#3 Cone Socket Front Surface Mirror
1	67-696/4*	#4 Cone Socket Front Surface Mirror
1	SVET-23E/18P	Griplite® S6 #23 Explorer Probe
1	SVET-5004	Griplite® S6 Jacquette J33/1Ys Scaler
1	SVET-767	Griplite® S6 Morse Scaler
1	SVET-5201	Griplite® S6 Columbia 13/14 Universal Curette
1	SVET-5202	Griplite® S6 Barnhart 1/2 Universal Curette
1	DTST**	Delron Sharpening Test Stick
1	STN-ARK-CON	Conical Arkansas Sharpening Stone
1	STN-ARK-FLT-4	#4 Flat Arkansas Sharpening Stone
1	3-084107	Storage Cassette 8" (20.3 cm) x 4-7/8" (10.1 cm) x 1-1/8" (2.5 cm)
1	-	Training DVD & Instructional Guide

*(not sold individually - reference is 12/bx) **(not sold individually - reference is 6/pk)



6806 - Extraction Instrument Kit

Kit includes:

Qty.	Reference	Descriptions
1	4-7	# 3 Scalpel Handle
1	73-64	Baby Molt Periosteal Elevator 6"
1	DEFF	Extracting Forceps Universal Mini 4-1/2"
1	5D-324	La Grange Scissors, Curved 4-1/4" One Serrated Blade Angled Shank
1	8-15TC	Olsen-Hegar Needle Holder 5-1/2" Tungsten Carbide
1	6-124	Brown-Adson Forceps 4-3/4" (7 x 7 Teeth)
1	5-16	Operating Scissors 5-1/2" Straight S/B
1	DELLUXS2	Luxating Elevator 4-1/2" Straight 2 mm
1	DELLUXS3	Luxating Elevator 4-1/2" Straight 3 mm
1	DELLUXS5	Luxating Elevator 4-1/2" Straight 5 mm
1	DELLUXC3	Luxating Elevator 4-1/2" Curved 3 mm
1	DELLUXC5	Luxating Elevator 4-1/2" Curved 5 mm
1	3-085009	Storage Cassette 8" (20.3 cm) x 6" (15.2 cm) x 1-3/16" (2.5 cm)



70-177 - Periotome Kit

Kit Includes:

Qty.	Reference	Description
1	70-170	Periotome Straight Double Ended
1	70-171	Periotome Hatchet Double Ended
1	70-172	Periotome Hoe Double Ended
1	3-072007	Storage Cassette 7" (17.7 cm) x 2-5/8" (5 cm) x 11/16" (1.9 cm)



* All buy/get promotions apply only to Integra™ Miltex® products advertised in this quarterly flyer unless otherwise noted.

* Redeemed goods must be of equal or lesser value of lowest priced purchased item.

Scissors (Surgical, Stitch and Bandage)

Buy 4, get 1 (mix and match)*

Description	Premium (TC)**	Premium (SS)**	Standard (SS)**
Mayo 5-1/2" Curved	5-122TC	5-122	MH5-122
Mayo 5-1/2" Straight	5-120TC	5-120	MH5-120
Mayo 6-3/4" Curved	5-126TC	5-126	MH5-126
Mayo 6-3/4" Straight	5-124TC	5-124	MH5-124
Metzenbaum 7" Curved	5-182TC	5-182	MH5-182
Metzenbaum 7" Straight	5-181TC	5-181	MH5-181
Metzenbaum 5" Curved (Delicate)	-	5-284	-
Metzenbaum 5" Straight (Delicate)	-	5-283	-
OR 5-1/2" Straight B/B	5-26TC	5-26	MH5-26
OR 5-1/2" Straight S/B	5-16TC	5-16	MH5-16
OR 5-1/2" Straight S/S	5-6TC	5-6	MH5-6
OR 5-1/2" Curved B/B	5-56TC	5-56	MH5-56
OR 5-1/2" Curved S/B	5-46TC	5-46	MH5-46
OR 5-1/2" Curved S/S	5-36TC	5-36	MH5-36

** TC = Tungsten Carbide

** SS = Stainless Steel



5-122TC

Needle Holders

Buy 4, get 1 (mix and match)*

Description	Ceram-A-Grip®	Premium (TC)**	Standard (SS)**
Olsen-Hegar with Suture Scissors 5-1/2" (Serrated)	8-15CG	8-15TC	MH8-15TC
Olsen-Hegar with Suture Scissors 6-1/2" (Serrated)	8-16CG	8-16TC	MH8-16TC
Olsen-Hegar with Suture Scissors 7-1/2" (Serrated)	-	8-17TC	MH8-17TC
Derf 4-3/4" (Serrated)	8-4CG	8-4TC	MH8-4TC
Webster 5-1/4" (Smooth)	8-6CG	8-6TC	MH8-6TC
Webster 4-3/4" (Smooth) (Extra Delicate)	8-7CG	8-7TC	-

** TC = Tungsten Carbide



8-15TC

Ring Handled Forceps

Buy 4, get 1 (mix and match)*

Description	Premium (SS)**	Standard (SS)**
Halsted Mosquito 5" Curved	7-4	MH7-4
Halsted Mosquito 5" Straight	7-2	MH7-2
Hartman Mosquito 3-1/2" Curved	7-26	MH7-26
Hartman Mosquito 3-1/2" Straight	7-24	MH7-24
Crile 5-1/2" Curved	7-44	MH7-44
Crile 5-1/2" Straight	7-42	MH7-42
Kelly 5-1/2" Curved	7-38	MH7-38
Kelly 5-1/2" Straight	7-36	MH7-36
Rochester-Carmalt 6-1/4" Curved	7-172	MH7-172

** SS = Stainless Steel



7-4

* All buy/get promotions apply only to Integra™ Miltex® products advertised in this quarterly flyer unless otherwise noted.

* Redeemed goods must be of equal or lesser value of lowest priced purchased item.

Dressing and Tissue Forceps Buy 4, get 1 (mix and match)*

Description	Premium (TC)**	Premium (SS)**	Standard (SS)**
Adson Dressing 4-3/4" (Serrated)	6-118TC	6-118	MH6-118
Adson Tissue 4-3/4" (1 x 2 Teeth)	-	6-120	MH6-120
Adson Tissue 4-3/4" (1 x 2 Teeth and Cross Serrated)	-	6-121	MH6-121
Adson Tissue & Suture 4-3/4" (1 x 2 Teeth with Tying Platform)	6-123TC	6-123	MH6-123
Brown-Adson Tissue 4-3/4" (7 x 7 Side Grasping Teeth)	6-124TC	6-124	MH6-124
Tissue 5" (1 x 2 Teeth)	-	6-42	MH6-42
Tissue 5-1/2" (1 x 2 Teeth)	-	6-44	MH6-44

** TC = Tungsten Carbide

** SS = Stainless Steel



6-118TC

Buy 2 General Surgery Kits, get \$100 in Sharp-It™ Cards and 1 box of Silitone™ Bands Free

6800 - General Surgery Kit

Kit includes:

Qty.	Reference	Description
1	5-124	Mayo Scissors 6-3/4" Straight
1	4-7	#3 Scalpel Handle
1	6008	Snook Hook 8"
3	7-4	Mosquito Forceps 5" Curved
2	7-36	Kelly Forceps 5-1/2" Straight
1	5-182	Metzenbaum Scissors 7" Curved
1	6-124	Brown-Adson Forceps 4-3/4" (7 x 7 Teeth)
4	7-172	Rochester-Carmalt Forceps 6-1/4" Curved
4	7-506	Backhaus Towel Clamps 5-1/4"
1	8-15TC	Olsen-Hegar Needle Holder 5-1/2" Tungsten Carbide
1	3-6815	Storage Cassette 15" (38.1 cm) x 8-1/4" (21 cm) x 1-1/4" (3.1 cm)

	Reference
Premium Grade	6800
Standard Grade	MH6800
Mid Grade	V96800



AHSHARPITCERT



VST

* All buy/get promotions apply only to Integra™ Miltex® products advertised in this quarterly flyer unless otherwise noted.

* Redeemed goods must be of equal or lesser value of lowest priced purchased item.

Mail or fax this completed coupon along with copy of dealer invoice dated between January 1 to June 30, 2012 to:
Attn: Coupon Redemption
Integra ▪ 589 Davies Drive, York, PA 17402 ▪ 800-229-6589 fax

Customer Service: 866-854-8300 USA
integralife.com/integra-miltex

List Your Redeemed Goods

Quantity	Reference #
_____	_____
_____	_____
_____	_____

Please send Integra™ Miltex® Products to:

Name of Person Ordering _____
Practice Name _____
Address _____
City/State/Zip _____
Phone _____
Dealer/Branch _____
Dealer Representative _____

Offer valid 1/1/12 - 6/30/12 in the USA only from end-users and does not apply to government, institutional or dealer stocking orders. Promotion void if copy of invoice is not attached with promotion coupon. Offers cannot be combined unless invoice totals meet promotion criteria. Coupon must be submitted within 30 days of expiration date of offer. Please allow up to 8 weeks for goods to be delivered. By participating in this offer, Customer agrees, as applicable, to properly disclose the discounts or reductions in price in the costs claimed or charges made under any federal or state health care program which provides cost or charge based reimbursement to Customer for the items covered under this offer.

- * Redeemed goods must be of equal or lesser value of lowest priced purchased item.
- * Redemption invoices subject to verification by Integra Miltex.
- * Buy/get promotion applies only to Integra™ Miltex® products advertised in this biannual flyer unless otherwise noted.



For more information or to place an order, please contact:
Integra ▪ 589 Davies Drive, York, PA 17402
866-854-8300 USA ▪ 717-840-2763 outside USA ▪ 717-840-9335 fax
integralife.com/integra-miltex

